



# Reserve Manifolds

Acme Reserve Manifolds connect reserve cylinders to deliver product as needed. They are available with stands or brackets for floor or wall mounting. When teamed with the Acme SLIQ, it provides continuous product delivery, as automatic switchover from the primary to secondary supply is achieved when primary dewars drop to a preset level.

## Features

- Available in single or doubled sided configurations
- Designs to accommodate 3 to 24 cylinders
- Includes individual cylinder supports per NFPA99 2002



## Ordering Information

### HP - Medical Reserve

| Gas Service             | CGA Inlet Fitting | Length | # of Cylinders | Part Number   |
|-------------------------|-------------------|--------|----------------|---------------|
| Oxygen                  | 540               | 48     | 10             | 272-701054-CV |
| Argon, Helium, Nitrogen | 580               | 48     | 10             | 272-781058-CV |

### Non-Medical Reserve

| Gas Service              | CGA Inlet Fitting | Length | # of Cylinders | Part Number |
|--------------------------|-------------------|--------|----------------|-------------|
| Oxygen                   | 540               | 48     | 10             | 272-701054  |
| Carbon Dioxide, Hydrogen | 350               | 48     | 10             | 272-761035  |
| Argon, Helium, Nitrogen  | 580               | 48     | 10             | 272-781058  |